



本PDF文件由 www.ichemistry.cn 免费提供, 全部信息请点击[215247-95-3](http://www.ichemistry.cn), 若要查询其它化学品请登录[CAS号查询网](http://www.ichemistry.cn)


如果您觉得本站对您的学习工作有帮助, 请与您的朋友一起分享:) [爱化学www.ichemistry.cn](http://www.ichemistry.cn)

CAS Number:215247-95-3 基本信息

中文名:	颜料紫20; 永固紫RL
英文名:	Diindolo[2,3-c:2',3'-n]triphenodioxazine, 9,19-dichloro-5,15-diethyl-5,15-dihydro-

C. I. Pigment Violet 23 (8CI);
 C. I. 51319;
 CFP-FF 802V;
 Carbazole Dioxazine Violet;
 Carbazole Violet;
 Carbazole Violet 23;
 Chromofine Violet 6510PK;
 ChromofineViolet RE;
 Cosmenyl Violet RL;
 Creanova 877-8895;
 Cromophtal Violet GT;
 Cyanadur Violet;
 Dioxazine Violet;
 Dioxazine purple;
 EB Violet 4B7906;
 EMCViolet RL 10;
 Fast Violet RL 3R Base;
 Fastogen Super Violet RN;
 Fastogen SuperViolet RN-S;
 Fastogen Super Violet RN-SU 02;
 Fastogen Super Violet RTS;
 Fastogen Super Violet RVS;
 Fastogen Super Violet RXE;
 Fastogen Super VioletRXS;
 Heliofast Red Violet EE;
 Heliofast Violet BN;
 Heliogen Violet;
 HeliogenViolet R Toner;
 Hostaperm RL-NF;
 Hostaperm Violet BL;
 Hostaperm Violet P-RL;
 Hostaperm Violet RL;
 Hostaperm Violet RL 02;
 Hostaperm Violet RL Special;
 Hostaperm Violet RL Special 14-4007;
 Hostaperm Violet RL Special TS;
 Hostaperm VioletRL-COF-VP 2777;
 Hostaperm Violet RL-NF;
 Hostaperm Violet RL-SP;
 HostapermViolet RL-SPL;

别名:
Lake Fast Violet RL;
Lake Fast Violet RLB;
Lionogen Violet HR;
Lionogen Violet R 6100;
Lionogen Violet R 6200;
Lionogen Violet RL;
LionolViolet HR;
Monolite Fast Violet R;
Novotex Violet BL VP 2435;
Novotex VioletBL-PC VP 2429;
Oriental Fast Violet BL;
PV 23;
PV Fast Violet BL;
PV FastViolet BLP;
PV Fast Violet RL;
PV Fast Violet RL-SPE;
Paliogen Violet 5890;
Paliogen Violet L 5890;
Permanent Violet;
Permanent Violet R;
Permanent VioletRL;
Pigment Violet 23;
Sandorin Violet BL;
Sanyo Permanent Violet BL-D 422;
Sumikacoat Fast Violet RSB;
Sumitomo Fast Violet RL 3R Base;
Sumitomo FastViolet RL Base;
Sumitone Fast Violet RL;
Sumitone Fast Violet RL 4R;
SumitoneFast Violet RL 4R Base;
Sumitone Fast Violet RL Base;
Sumitone Fast Violet RLS;
Sun Violet VHD 6003;
Sunfast Violet 23;
Symuler Fast Violet BBL;
Symuler FastViolet BBLN;
Unisperse Violet B-E;
Unisperse Violet B-S;
VHD 6003;
Violet P-RL;
Viscofil Violet BLN;
Vynamon Violet 2B

分子结构:	
分子式:	C ₃₄ H ₂₂ Cl ₂ N ₄ O ₂
CAS登录号:	215247-95-3
物理化学性质	
性质描述:	颜料紫20 (215247-95-3) 的性状: 蓝光紫色粉末, 色泽鲜艳。具有高的着色力, 优异的耐热性、耐渗化性和耐晒性。 详情请看
CAS#215247-95-3化学试剂供应商(点击生产商链接可查看价格)	
<p style="color: red;">供应商信息已更新, 请登录爱化学 CAS No. 215247-95-3 查看</p> <p>若您在此化学品供应商, 请按照化工产品收录说明进行免费添加</p>	
其他信息	
产品应用:	颜料紫20 (215247-95-3) 的用途: 主要用于涂料、油墨、橡胶和塑料的着色, 也用于合成纤维的原液着色。安全性10kg纸板桶或铁桶包装, 内衬塑料袋。 详情请看
生产方法及其他:	颜料紫20 (215247-95-3) 的制备方法: 3-氨基-9-乙基咪唑、邻二氯苯、水、醋酸钠、苄基三甲氯化铵与四氯苯醌的混合物在40~60℃加热1h。用苯磺酰氯处理并加热至165~170℃, 在此温度维持6h, 产品收率可达84%。然后再升温至150~160℃, 蒸出反应物中的水。 详情请看
相关化学品信息	
216144-77-3 218773-48-9 211991-79-6 211487-66-0 2194-21-0 2199-64-6 3-羧基吡嗪 21879-83-4 21948-10-7 2165-21-1 219959-85-0 21121-54-0 21250-85-1 反式-α-氰基肉桂酸乙酯 215777-46-1 437	
生成时间2015-2-2 12:26:52	